

USE:

1. Include skin inspections, proper nutrition, turning schedules and daily inflation checks as part of a balanced wound care management program.
2. This product is intended for use as a medical device to aid in the prevention and treatment of pressure ulcers or “bed sores” and general comfort.
3. Weight limit for 200WC/218WC/228WC: Up to 300lbs. depending on body type.
4. Weight limit for 240WCI: Up to 700lbs. depending on body type.
5. Weight limit for 214WCI: Up to 100lbs. depending on body type.

MAINTENANCE:

1. Products should be at normal room temperature before inflating.
Store in room temperature of 55 to 85 degrees.
2. Clean all WAFFLE products with soapy water or bleach/water. Do not use cleaning solvents containing alcohol or petroleum distillates. Refer to infectious disease protocol and follow disinfectant manufacturers’ instructions. Do not wash or dry in institutional laundry. Proper product inflation and skin care protocol must be maintained to maximize the skin care management provided by these products. Skin inspection, daily hand checks and turning schedules are vital to your skin care program.
3. Use facility procedure for disposing of product and packaging materials.
4. Altitude and temperature may affect inflation.

PLEASE READ:

1. An anti-microbial agent is blended within the material formula of the product.
2. This product has flame resistant qualities in accordance with California Technical Bulletin 129 and ASTM E1590.
3. Do not smoke or use near open flame.
4. This product is not sterile and is not intended to be sterilized.
5. Keep packaging materials away from children.
6. All WAFFLE products are “**Natural Rubber Latex**” FREE.
7. All EHOB products are backed by a limited warranty that varies in duration depending on the product. Products found defective in workmanship or materials will be replaced provided failure is not the result of abuse or misuse, alternation or modification, and occurs within the specified warranty period. Tampering with product voids warranty; Example: glue or any unauthorized patching device. Please contact our customer service department for full details. Product warranties are in lieu of warranties for merchantability fitness for a particular purpose and all other warranties expressed or implied.



The Clinical Leader of Medical Devices for Patient Care

250 North Belmont Avenue, Indianapolis, IN 46222
Ph. **800.899.5553**
Fx. **317.972.4625**

Web: www.ehob.com, E-Mail: corporate@ehob.com
ISO 9001: 2008



WAFFLE® Brand Cushions



Indications for use:

Comfort, prevention and healing therapy through Stage IV pressure ulcers. Deep tissue injury protection.

Product	Model	Size	Warranty	WAFFLE L.A.D. Pump Inflation
WAFFLE Original Cushion	200WC 200WCI	17" x 17" x 1.5"	90 Day	9 Complete Strokes
WAFFLE Standard Adult Cushion	218WC 218WCI	19" x 19" x 1.5"	90 Day	14 Complete Strokes
WAFFLE TPU Cushion	228WC 228WCI	19" x 19" x 1.5"	6 Month	14 Complete Strokes
WAFFLE Bariatric Cushion	240WCI	22" x 28" x 2"	180 Day	Shipped Inflated
WAFFLE Pediatric Cushion	214PCI	14" x 14"	90 Day	Shipped Inflated

Cushions are shipped pre-inflated.

Cushion covers available for the 200WC, 218WC and 240WCI Single Patient Use



WAFFLE Cushion Hand Check Procedure



Keep your hand flat, do not bend fingers under the WAFFLE Cushion.

While patient is seated on the WAFFLE Cushion, slide your hand, palm side up, under the cushion and between patient's legs.

The ability to lift the fingers upward before coming into contact with the patient's body is an indicator of having the correct amount of air.

You may feel the patient as your hand has displaced the air that has been supporting their weight.

The purpose of the hand check is to ensure there is enough air in the cushion, thereby providing pressure redistribution for the patient.

Remember to perform hand checks every shift and to encourage the patient to shift their weight every 15 minutes.



Most WAFFLE Cushions are shipped inflated. However, **IF** you are using a hand pump to inflate your product, follow the inflation instructions listed below.

Unfold the product to expose inflation valve.


Using the WAFFLE L.A.D Pump insert the tip of the pump into the valve.



WAFFLE L.A.D. Hand Pump

Refer to the inflation chart on the side of the pump and inflate your model using the recommended level of strokes.

Important: One out and in stroke equals One Stroke.



Cushion is properly inflated when you can easily roll one side just past the first set of holes.

Six Quick Tips on Your Cushion

1. The WAFFLE Cushion will appear only half full to achieve its full benefit.
2. Patient should sit **IN** the cushion and not **ON** it.
3. You should be able to roll one side of the cushion past the first set of holes.
4. Placing the cushion in a pillow case provides extra coolness and comfort.
5. Cushion should be placed in chair with valve down and to the back of surface.
6. Reinflating your cushion may require less air.

Other Uses for Your WAFFLE Cushion

Between Knees



Protect Ulnar Nerve



Protect Occiput

